

20Cb-3 Stainless

Corrosion Resistance of 20Cb-3 Stainless

The corrosion resistance of 20Cb-3 stainless in pure sulfuric acid at 176°F (80°C) and at the boiling point are shown in Figures 1 and 2. The corrosion resistance is expressed as a corrosion rate in inches penetration per year (ipy). The corrosion rates in pure sulfuric acid should be used only as a guide since many contaminants in commercial sulfuric acid decrease the corrosive attack. Corrosion rates in sulfuric acid are actually much lower in the presence of iron, copper, and chromium ions, usually present in pickling and plating solutions. Higher rates may occur in the presence of chloride ions.

At 176°F (80°C) 20Cb-3 stainless exhibits a slight improvement over the former lower nickel grade, Carpenter 20 stainless. In the region of the "hump," which exists in the range of 60 to 75% sulfuric acid, 20Cb-3 stainless generally exhibits about half the corrosion rate experienced in the lower nickel grade. Although 20Cb-3 stainless would not usually be recommended for service in 60 to 75% sulfuric acid, the alloy offers an advantage in the case of accidental dilution of acid into this range.

The corrosion rates in boiling sulfuric acid for 20Cb-3 stainless are significantly lower than the lower nickel grade -- particularly above 10% acid.

Figure 1




Figure 2




Figure 3 shows an iso-corrosion chart for 20Cb-3 stainless. The chart is divided into zones according to the corrosion resistance of 20Cb-3 stainless exposed to the range of temperatures and concentrations represented by each zone. These data are based on laboratory tests in reagent grade H_2SO_4 with no intentional aeration or deaeration. The chart provides resistance to only general corrosion and is intended for general guidance as the corrosion resistance can vary with impurities, aeration, etc. See page 13 for other factors that influence corrosion resistance.

Heat Transfer Conditions

The corrosion rates presented in Figure 4 were obtained under conditions that are frequently encountered in the chemical process industries. Many times the application of a material involves a situation where the metal is maintained at a temperature above that of the corrosion environment. Thus, heat is transferred from the metal to the solution, with the metal acting as a heat transfer surface. Corrosion rates under these conditions can be markedly higher than those experienced where the metal and solution are at the same temperature and atmospheric pressure. These heat transfer tests were conducted by sealing a circular specimen to the bottom of a custom designed glass vessel.

Figure 3





Figure 4



20Cb-3 Stainless

One side of the specimen was exposed to the corrodent while the other side was heated by an external source. This enabled testing to be conducted where the metal specimen temperature was higher than the temperature of the corrodent, simulating heat exchanger conditions.

The curves shown in Figure 4 indicate the corrosion rate of 20Cb-3 stainless to various sulfuric acid concentrations as a function of the metal temperature. While the sulfuric acid solutions were maintained at their normal atmospheric boiling points of 214/223°F (101/106°C), corrosion rates were obtained at metal temperatures of 221°F (105°C) to 257°F (125°C). These data are particularly useful where sulfuric acid solutions are being heated, such as in pickling tanks.

Stress-Corrosion Resistance

The nickel content of 20Cb-3 stainless has provided resistance to stress-corrosion cracking in sulfuric acid at a variety of temperatures and concentrations. 20Cb-3 stainless also has better resistance than lower-nickel austenitic stainless to chloride stress-corrosion cracking in magnesium chloride (MgCl_2) boiling at 311°F (155°C), a standard test for susceptibility to this cracking mechanism. While the alloy has been made to crack in this unusually severe test, it has provided excellent resistance in most service environments. For example, 20Cb-3 stainless did not crack in boiling acidified 25% sodium chloride (NaCl , pH 1.5), an environment which cracked Type 304 and Type 316 stainless steels. This sodium chloride test is believed to provide a better indication of service performance than boiling magnesium chloride solutions.

Intergranular Corrosion Resistance

In the annealed and annealed plus sensitized [1250°F (677°C), one hour] conditions, 20Cb-3 stainless will pass the nitric acid test (240 hours in boiling 65% nitric acid) and the ferric sulfate-sulfuric acid test (120 hours in a 50% sulfuric acid solution containing ferric sulfate) with a rate of 0.002 ipm or less.

The nitric acid test is described in ASTM A-262, Practice C. The ferric sulfate-sulfuric acid test is described in ASTM A-262, Practice B, and ASTM G-28, Method A.

The 24-hour copper accelerated acidified copper sulfate test in A-262, Practice E, is also applicable to 20Cb-3 stainless in the annealed and annealed plus sensitized [1250°F (677°C), one hour] condition.

See page 30 for the heat treatment to properly stabilize 20Cb-3 stainless.

Sour Service

The excellent resistance of 20Cb-3 stainless to sulfide stress cracking has enabled inclusion of this alloy in NACE MR0175, "Sulfide Stress Cracking Resistant Metallic Materials for Oil Field Equipment." The alloy (UNS N08020) is acceptable in the annealed or cold worked condition at a hardness level of Rockwell C 32 maximum.

Scaling

20Cb-3 stainless has a safe scaling temperature in continuous service of 1800°F (982°C). Figures 5 and 6 compare the oxidation resistance of 20Cb-3 stainless with that of Types 304 and 316 stainless steels. The data shown in these figures were developed in Carpenter's Research Laboratory using an electric furnace with a still air atmosphere. Oxidation rates in field service may vary depending on the atmosphere, temperature, etc.

Figure 5

Comparative Oxidation at Various Temperatures for 16 Hours





Figure 6

Comparative Oxidation at 1800°F (982°C) for Extended Periods of Time



20Cb-3 Stainless

Corrosion Resistance Table

This table shows the resistance of a number of materials to the more common chemicals. Many factors influence the resistance of materials to various solutions. Factors which must be given consideration for service in corrosive environments are: temperature, concentration, aeration, influence of inhibiting or accelerating contaminants, influence of recirculation, solids in suspension, velocity, frequency of use, and equipment design. See page 13 for details.


The influence of contaminants is probably the most important from a commercial standpoint. Few corrosive solutions will be free of all contaminants. The majority of these contaminants have no influence on corrosion, but those that do generally affect the conditions greatly.

The corrosion data for all grades except 20Cb-3 stainless is reprinted from Corrosion Data Survey, 1974 Edition, published by the National Association of Corrosive Engineers. The corrosion rates for 20Cb-3 stainless represent a composite of the NACE Corrosion Data Survey and more current data developed in Carpenter's Corrosion laboratory.

Code:

- Corrosion Rate less than 0.002" per year
- Corrosion Rate less than 0.020" per year
- Corrosion Rate from 0.020" to 0.050" per year
- X Corrosion Rate greater than 0.050" per year

Because the in-service performance of stainless steels may be significantly affected by minor variations in environment and application, Carpenter Specialty Alloys makes no warranty, either express or implied, and accepts no liability, compensatory or consequential, for the performance of stainless steels in specific applications for which this publication may be the basis for information. See disclaimer on inside front cover.



Corrodent	Carpenter 20Cb-3 Stainless	Type 316	Hastelloy* Alloy B	Hastelloy* Alloy C	Brass Cu 70/80%	Monel**
Acetic Acid Aerated 7, 12						
Acetic Acid No Air						
Acetic Acid Vapor						
Acetic Anhydride In Acetic Acid 4, 7, 9						
Acetone 4						
Aluminum Chloride 8, 13						
Aluminum Potassium Sulfate						
Aluminum Sulfate						
Ammonium Carbonate						
Ammonium Chloride						
Ammonium Nitrate 3, 4						
Ammonium Sulfate						
Amyl Acetate 4						
Aniline 1, 7, 9						
Aniline Sulfite						
Arsenic Acid 1						
Barium Carbonate 1						

Corrodent	Carpenter 20Cb-3 Stainless	Type 316	Hastelloy* Alloy B	Hastelloy* Alloy C	Brass Cu 70/80%	Monel**
Barium Chloride 1						
Barium Hydroxide 1						
Barium Nitrate 1, 4						
Barium Sulfate						
Benzene 1, 2, 3, 4, 14						
Benzene Sulfonic Acid 2						
Benzoic Acid						
Boric Acid						
Butyl Acetate 4, 7, 9						
Butyric Acid						
Cadmium Sulfate						
Calcium Bisulfite						
Calcium Carbonate						
Calcium Chlorate 4						
Calcium Chloride						
Calcium Sulfate						
Camphor						

Corrodant	Carpenter 20Cb-3 Stainless	Type 316	Hastelloy* Alloy B	Hastelloy* Alloy C	Brass Cu 70/30%	Monel**
Carbonic Acid 6, 9	20 [Corrosion pattern]	50C [Corrosion pattern]	[Corrosion pattern]	[Corrosion pattern]	4 2 [Corrosion pattern]	40 41 [Corrosion pattern]
Carbon Disulfide 1, 4, 9, 15	[Corrosion pattern]	[Corrosion pattern]	[Corrosion pattern]	[Corrosion pattern]	[Corrosion pattern]	[Corrosion pattern]
Carbon Monoxide 1, 4	60C [Corrosion pattern]	50C [Corrosion pattern]	100C [Corrosion pattern]	40C [Corrosion pattern]	50C 17 [Corrosion pattern]	42 33 [Corrosion pattern]
Carbon Tetra- chloride 1, 9	1 [Corrosion pattern]	5 1 [Corrosion pattern]	[Corrosion pattern]	[Corrosion pattern]	1 [Corrosion pattern]	[Corrosion pattern]
Chlorine 1, 16	[Corrosion pattern]	4 5 [Corrosion pattern]	[Corrosion pattern]	[Corrosion pattern]	1 2 43 [Corrosion pattern]	[Corrosion pattern]
Chloroform 1, 9	[Corrosion pattern]	4 5 1 [Corrosion pattern]	[Corrosion pattern]	[Corrosion pattern]	[Corrosion pattern]	[Corrosion pattern]
Chloro- sulfonic Acid 1, 7, 11	[Corrosion pattern]	4 [Corrosion pattern]	[Corrosion pattern]	[Corrosion pattern]	[Corrosion pattern]	[Corrosion pattern]
Chromic Acid 7, 10	44 [Corrosion pattern]	8 44 [Corrosion pattern]	[Corrosion pattern]	[Corrosion pattern]	[Corrosion pattern]	5 [Corrosion pattern]
Chromic Sulfates	[Corrosion pattern]	[Corrosion pattern]	[Corrosion pattern]	[Corrosion pattern]	2 [Corrosion pattern]	[Corrosion pattern]
Citric Acid	[Corrosion pattern]	4 45 [Corrosion pattern]	[Corrosion pattern]	[Corrosion pattern]	[Corrosion pattern]	[Corrosion pattern]
Copper Nitrate 1, 10	[Corrosion pattern]	[Corrosion pattern]	[Corrosion pattern]	[Corrosion pattern]	[Corrosion pattern]	[Corrosion pattern]
Copper Sulfate 1	[Corrosion pattern]	9 10 [Corrosion pattern]	[Corrosion pattern]	[Corrosion pattern]	2 47 [Corrosion pattern]	2 [Corrosion pattern]
Cupric Cyanide 1	[Corrosion pattern]	[Corrosion pattern]	[Corrosion pattern]	[Corrosion pattern]	[Corrosion pattern]	[Corrosion pattern]
Ethanol	[Corrosion pattern]	[Corrosion pattern]	31 [Corrosion pattern]	31 [Corrosion pattern]	[Corrosion pattern]	3 [Corrosion pattern]
Ethyl Acetate 4	[Corrosion pattern]	[Corrosion pattern]	[Corrosion pattern]	[Corrosion pattern]	2 [Corrosion pattern]	[Corrosion pattern]
Ethyl- Chloride 4	[Corrosion pattern]	1 5 [Corrosion pattern]	[Corrosion pattern]	[Corrosion pattern]	1 [Corrosion pattern]	[Corrosion pattern]
Ethylene Dibrom- ide 1, 8, 9	4 [Corrosion pattern]	4 [Corrosion pattern]	[Corrosion pattern]	[Corrosion pattern]	[Corrosion pattern]	1 [Corrosion pattern]

Corrodant	Carpenter 20Cb-3 Stainless	Type 316	Hastelloy* Alloy B	Hastelloy* Alloy C	Brass Cu 70/30%	Monel**
Ethylene Dichloride 4, 9	48 [Corrosion pattern]	21 1 [Corrosion pattern]	[Corrosion pattern]	[Corrosion pattern]	[Corrosion pattern]	[Corrosion pattern]
Ferric Nitrate 4	[Corrosion pattern]	[Corrosion pattern]	[Corrosion pattern]	[Corrosion pattern]	[Corrosion pattern]	[Corrosion pattern]
Ferric Sulfate	[Corrosion pattern]	[Corrosion pattern]	[Corrosion pattern]	[Corrosion pattern]	[Corrosion pattern]	[Corrosion pattern]
Ferrous Sulfate	[Corrosion pattern]	8 [Corrosion pattern]	[Corrosion pattern]	[Corrosion pattern]	2 [Corrosion pattern]	49 [Corrosion pattern]
Fluorine Gas 1, 4, 17	[Corrosion pattern]	[Corrosion pattern]	[Corrosion pattern]	[Corrosion pattern]	[Corrosion pattern]	250C 32 [Corrosion pattern]
Fluosilicic Acid 1	4 [Corrosion pattern]	4 6 [Corrosion pattern]	5 [Corrosion pattern]	5 [Corrosion pattern]	2 [Corrosion pattern]	5 50 [Corrosion pattern]
Formal- dehyde 2, 4, 6, 8	[Corrosion pattern]	20 [Corrosion pattern]	[Corrosion pattern]	[Corrosion pattern]	7 [Corrosion pattern]	7 [Corrosion pattern]
Formic Acid 3, 7, 9	[Corrosion pattern]	4 51 [Corrosion pattern]	52 [Corrosion pattern]	[Corrosion pattern]	2 [Corrosion pattern]	2 [Corrosion pattern]
Furfural 4	[Corrosion pattern]	[Corrosion pattern]	[Corrosion pattern]	[Corrosion pattern]	[Corrosion pattern]	[Corrosion pattern]
Gallic Acid	[Corrosion pattern]	[Corrosion pattern]	[Corrosion pattern]	[Corrosion pattern]	[Corrosion pattern]	[Corrosion pattern]
Glutamic Acid	[Corrosion pattern]	4 5 8 [Corrosion pattern]	54 [Corrosion pattern]	54 [Corrosion pattern]	4 7 [Corrosion pattern]	[Corrosion pattern]
Glycerol	4 53 8 [Corrosion pattern]	[Corrosion pattern]	[Corrosion pattern]	[Corrosion pattern]	[Corrosion pattern]	[Corrosion pattern]
Hexa- methylene Tetramine 4, 7	[Corrosion pattern]	[Corrosion pattern]	[Corrosion pattern]	[Corrosion pattern]	7 [Corrosion pattern]	[Corrosion pattern]
Hydro- chloric (Aerated) 1, 8, 7, 9	4 [Corrosion pattern]	[Corrosion pattern]	55 [Corrosion pattern]	[Corrosion pattern]	[Corrosion pattern]	10 [Corrosion pattern]
Hydro- chloric Acid (No Air)	4 [Corrosion pattern]	[Corrosion pattern]	55 [Corrosion pattern]	[Corrosion pattern]	[Corrosion pattern]	[Corrosion pattern]
Hydrocy- anic Acid + Hydrogen Cyanide 3, 4, 7, 9	14 [Corrosion pattern]	6 [Corrosion pattern]	56 57 [Corrosion pattern]	[Corrosion pattern]	7-17 [Corrosion pattern]	56 57 [Corrosion pattern]
Hydrofluoric Acid (Aerated) 1, 7, 9	[Corrosion pattern]	[Corrosion pattern]	[Corrosion pattern]	[Corrosion pattern]	[Corrosion pattern]	5 [Corrosion pattern]

Corrodent	Carpenter 20Cb-3 Stainless	Type 316	Hastelloy* Alloy B	Hastelloy* Alloy C	Brass Cu 70/80%	Monel**
Hydrogen	58 59 60 61 4	58 59 60 61 4	58 59 60 61 4	58 59 60 61 4	58 59 60 61 4	58 59 60 61 4
Hydrogen Chloride (Anhydrous)	7, 9	7, 9	7, 9	7, 9	7, 9	7, 9
Hydrogen Fluoride (Anhydrous)	63 64	63 64	63 64	63 64	63 64	63 64
Hydrogen Peroxide	3, 7	3, 7	3, 7	3, 7	3, 7	3, 7
Hydrogen Sulfide - Dry	1, 4	1, 4	1, 4	1, 4	1, 4	1, 4
Lactic Acid	4, 16 8 9	4, 16 8 9	4, 16 8 9	4, 16 8 9	4, 16 8 9	4, 16 8 9
Lead Acetate	1, 6	1, 6	1, 6	1, 6	1, 6	1, 6
Lead Nitrate	1, 10	1, 10	1, 10	1, 10	1, 10	1, 10
Lithium Chloride	1	1	1	1	1	1
Lithium Hydroxide	2	2	2	2	2	2
Magnesium Chloride	5	5	5	5	5	5
Magnesium Chloride + Calcium Chloride	17 20 4 6	17 20 4 6	17 20 4 6	17 20 4 6	17 20 4 6	17 20 4 6
Magnesium Hydroxide Or Magnesium Oxide	9	9	9	9	9	9
Magnesium Sulfate	10	10	10	10	10	10
Maleic Acid	1, 7	1, 7	1, 7	1, 7	1, 7	1, 7
Malic Acid	2	2	2	2	2	2
Manganese Chloride	5	5	5	5	5	5

Corrodent	Carpenter 20Cb-3 Stainless	Type 316	Hastelloy* Alloy B	Hastelloy* Alloy C	Brass Cu 70/80%	Monel**
Mercuric Nitrate	1	1	1	1	1	1
Mercury	5	5	5	5	5	5
Methane	3, 4	3, 4	3, 4	3, 4	3, 4	3, 4
Methanol	1, 4, 9	1, 4, 9	1, 4, 9	1, 4, 9	1, 4, 9	1, 4, 9
Methyl Chloride	1, 4, 9	1, 4, 9	1, 4, 9	1, 4, 9	1, 4, 9	1, 4, 9
Methylene Chloride	6, 7	6, 7	6, 7	6, 7	6, 7	6, 7
Mixed Acids $H_2SO_4 + HNO_3$	68	68	68	68	68	68
Mono- ethanol- amine	16	16	16	16	16	16
Naph- thalene	4	4	4	4	4	4
Naphthenic Acid	2 33	2 33	2 33	2 33	2 33	2 33
Nickel Chloride	2	2	2	2	2	2
Nickel Nitrate	4	4	4	4	4	4
Nickel Sulfate	2	2	2	2	2	2
Nitric Acid	1, 6, 7, 18 9 69 37	1, 6, 7, 18 9 69 37	1, 6, 7, 18 9 69 37	1, 6, 7, 18 9 69 37	1, 6, 7, 18 9 69 37	1, 6, 7, 18 9 69 37
Nitric Acid - Red Fuming	7, 9, 18	7, 9, 18	7, 9, 18	7, 9, 18	7, 9, 18	7, 9, 18
Nitric + Hydro- fluoric Acids						
Nitriding Gases						

Corrodent	Carpenter 20Cb-3 Stainless	Type 316	Hastelloy* Alloy B	Hastelloy* Alloy C	Brass Cu 70/80%	Monel**
Nitro- Benzene 1, 7, 9	22					
Nitro- glycerin 3						
Nitrous Acids						
Nitrous Oxide						
Oleic Acid 20					4, 7, 32	
Oxalic Acid 1, 7	9				2	3
Per- chloric Acid 7, 18, 19						
Per- chloro- ethylene 9						
Phenol 1, 7					7, 17	
Phosphoric Acid (Aerated) 7						
Phosphorus 1, 4, 7, 20						
Phthalic Anhydride 7						
Phthalic Anhydride (Pure) + Maleic Anhydride						
Picric Acid 1, 3, 4						
Potassium Bromide 13					2	2
Potassium Carbonate					7-17	
Potassium Chlorate 1, 10					2	

Corrodent	Carpenter 20Cb-3 Stainless	Type 316	Hastelloy* Alloy B	Hastelloy* Alloy C	Brass Cu 70/80%	Monel**
Potassium Chloride 4	4	4			5	
Potassium Chromate 5, 22						
Potassium Cyanide 1, 6, 7			23		2, 4	2, 23
Potassium Dichromate 1, 7						
Potassium Ferri- cyanide 1		71				
Potassium Ferro- cyanide 1						
Potassium Hydroxide 1, 5, 7, 8, 21	2, 19	5			2, 7	2, 5
Potassium Nitrate 4, 10						
Potassium Oxalate						
Potassium Permanganate 4, 10						
Potassium Peroxide 3, 4, 10						
Potassium Sulfate						
Propionic Acid 2						2
Pyridine 4, 9					2, 7	
Pyrogalllic Acid						
Pyro- lignous Acid 2						
Rosin 4						

Corrodent	Carpenter 20Cb-3 Stainless	Type 316	Hastelloy* Alloy B	Hastelloy* Alloy C	Brass Cu 70/80%	Monel**
Salicylic Acid 22						
Silver Nitrate 1.7						
Sodium Acetate						
Sodium Aluminum Sulfate						
Sodium Bicar- bonate						
Sodium Bichromate 5, 6						
Sodium Bisulfate						
Sodium Bisulfite						
Sodium Bromide						
Sodium Carbonate						
Sodium Chlorate 4						
Sodium Chloride						
Sodium Chromate						
Sodium Citrate						
Sodium Cyanide 1						
Sodium Ferri- cyanide 1						
Sodium Formal- dehyde Sulfoxylate						

Corrodent	Carpenter 20Cb-3 Stainless	Type 316	Hastelloy* Alloy B	Hastelloy* Alloy C	Brass Cu 70/80%	Monel**
Sodium Hydro- sulfide 1.7.9						
Sodium Hydroxide 1.6						
Sodium Meta- silicates						
Sodium Nitrate 10						
Sodium Perborate 10						
Sodium Perchlorate 75						
Sodium Phosphate 10 14						
Sodium Phosphate (Tribasic) 5 - 17						
Sodium Silicates 600 - 1000						
Sodium Sulfate 3						
Sodium Sulfide 9						
Sodium Sulfite						
Stannous Chloride 4						
Stearic Acid						
Sulfate Black Liquor						
Sulfate Green Liquor						
Sulfite Liquor with 10% Sulfur Dioxide 4						

Footnotes

Footnotes for Corrosives:

- 1** Poison
- 2** Toxic
- 3** Explosive
- 4** Flammable
- 5** Ingestion poison
- 6** Inhalant poison
- 7** Attacks skin
- 8** Irritant
- 9** Vapor harmful
- 10** Ignites organics
- 11** Fuming liquid
- 12** Hygroscopic
- 13** Liberates HCl in water
- 14** Narcotic
- 15** Volatile
- 16** Hazardous under pressure
- 17** Ignites combustibles
- 18** Fire hazard
- 19** Explosive over 70%
- 20** Ignites in moist air at 30°C
- 21** Exothermic in water
- 22** Dust explodes
- 23** Explosive dust
- 24** Exothermic with water

Footnotes for Data Squares:

- 1** No water
- 2** No air, oxygen
- 3** Low air, oxygen
- 4** Pits
- 5** Stress cracks
- 6** Stress corrosion
- 7** Discolors
- 8** Crevice attack
- 9** Intergranular attack
- 10** No chlorides
- 11** May discolor
- 12** May catalyze
- 13** May pit
- 14** May stress crack
- 15** Transgranular attack
- 16** Vapor
- 17** Aerated
- 18** Catalyzes
- 19** Static
- 20** Agitated
- 21** ~7 pH
- 22** <7 pH
- 23** >7 pH
- 24** No HCl, H₂SO₄, NaCl
- 25** No ferric chloride
- 26** ~0.1% acetic acid
- 27** Also sludge
- 28** No iron salts
- 29** No sulfuric acid
- 30** Explosive
- 31** With H₂SO₄
- 32** With steam
- 33** No sulfur
- 34** No stress
- 35** No ammonia
- 36** 300 psi
- 37** Stress relieved
- 38** No HCl, Cu, Ni ions
- 39** No Cu, Fe ions
- 40** Over 70% air
- 41** 20-70% air, 530 psi
- 42** With sulfur, <340°C = x
- 43** <10 mg/l
- 44** No H₂SO₄
- 45** <60 psi
- 46** No sulfides
- 47** <20% zinc
- 48** Trace HCl
- 49** pH 2 to 3.5
- 50** Annealed, immersed
- 51** >2.25% Mo
- 52** Erratic
- 53** With NaCl
- 54** With NaCl, HCl, H₂O₂
- 55** No Fe, Cl
- 56** With + ~0.05-1% H₃PO₄ or H₂SO₄
- 57** +SO₂ or HCOOH
- 58** >RC 22, 60,000
- 59** Annealed
- 60** No cold work
- 61** No H₂S
- 62** Permeable to H₂
- 63** Unsulfated
- 64** With or without steam
- 65** 240 psi
- 66** Cold worked
- 67** >80% copper
- 68** >20% sulfuric, bal. nitric acid
- 69** No Mo; low C
- 70** Red fuming
- 71** Pits in chlorides
- 72** Over 400°C
- 73** Steam and air
- 74** 75-100% concentration
- 75** Low NaCl
- 76** With HCl
- 77** <17% zinc
- 78** <0.23%, 200 psi
- 79** 300 psi
- 80** No SO₃
- 81** No SO₃
- 82** High pressure
- 83** 75-120 psi
- 84** No sodium sulfite
- 85** + ammonia
- 86** Avoid hydroxides
- 87** Saturated

* Registered Trademark of
Haynes International, Inc.

** Registered Trademark of INCO