

## Appendix 2 Chemical Compatibility of Materials with Homogeneous GRAFOIL Sheet

Graphite is one of the most thermally stable and chemically resistant materials known to man. GRAFOIL flexible graphite is manufactured employing a process that introduces no organic or inorganic elastomeric binders, fillers, or other potentially fugitive ingredients that could limit its thermal and chemical stability.

GRAFOIL flexible graphite is resistant to attack from nearly all organic and inorganic fluids with the exception of highly oxidizing chemicals and concentrated, highly oxidizing mineral acids such as nitric or sulfuric acids. GRAFOIL flexible graphite is an extremely chemically resistant material. The classes of organic chemicals that should not be used with GRAFOIL flexible graphite are those that are highly oxidizing, such as nitrates, persulfates, perbenzoate, and peroxides. Unacceptable compatibility for inorganic chemicals would include molten sodium or potassium hydroxide and chlorine dioxide.

If these chemicals are used in a critical or dangerous process, we always recommend their evaluation in a test loop since we do not have control over the environment or process in which the GRAFOIL flexible graphite gaskets are used.

It is much easier to list the chemicals which GRAFOIL flexible graphite is not compatible with instead of listing all of the chemicals it is compatible. Table XVI lists the strong oxidizers in which GRAFOIL flexible graphite would not be compatible and Table XVII lists the vast majority of chemicals which are compatible and could possibly come in contact with GRAFOIL flexible graphite during industrial use.

### Table XVI Strong Oxidizers

(not recommended with GRAFOIL flexible graphite)

Aqua Regia	Bromine (dry)	Calcium Chlorate
Calcium Hypochlorite*	Calcium Nitrate	Chlorazotic Acid
Chlorine Dioxide*	Chlorine Trifluoride	Chloric Acid
Chloroazotic Acid	Chloronitrous Acid	Chromates
Chromic Acid	Chromic Anhydride	Chromium
Chromium Trioxide*	Dichloropropionic Acid*	Dichromates
Hydrogen Dioxide	Hydrogen Peroxide**	Lime Nitrate
Lime Saltpeter	Molten Alkaline	Nitrates
Nitric Acid*	Nitric Oxide	Nitrocalcite

Nitrohydrochloric Acid	Nitromuriatic Acid	Norge Niter
Norwegian Saltpeter	Oleum (Fuming Sulfuric Acid)	Oxygen*(above +600°F)
Ozone	Perchloric Acid	Permanganate solutions
Persulfates	Perbenzoates	Perborates
Peroxide	Potassium Bichromate	Potassium Chlorate
Potassium Chromate	Potassium Dichromate	Potassium Nitrate
Sodium Chlorite*(over 4%)	Sodium Hypochlorite*	Sodium Peroxide
Sulfuric Acid	Sulfur Trioxide*	

\*GRAFOIL flexible graphite may be acceptable in specific temperature/concentration ranges. Call Graftech to discuss your specific application.

\*\* See special limitation notes on page 30.

Table XVII is a listing for many of the chemicals which are chemically compatible with GRAFOIL flexible graphite. The information has been taken from various sources of information.

Again, if these chemicals are used in a critical or dangerous process, we always recommend their evaluation in a test loop, since we do not have control over the environment or process in which the GRAFOIL gaskets are used.

Care must be taken of what combinations of GRAFOIL flexible graphite, interlayer material, chemicals, temperature and concentration are used together. The combination may not be compatible.

**Table XVII Chemicals which are compatible with GRAFOIL Flexible Graphite**

**Chemical Reagents**

Acetaldehyde	Alcohol
Acetate Solvents	Alcohol (Ethyl)
Acetic Acid	Aldehyde
Acetic Anhydride	Alkanes
Acetone	Alkyl Acetone
Acetone, Dry	Alkyl Alcohol
Acetylene Gas	Alkyl Amine
Acidaldehyde	Alkyl Arylsulphonics
Acrylonitrile	Alkyl Benzene
Air, Instrument	Alkyl Chloride
Air, Plant	Alkylate, Light
Air, Process	Alum
Air, Starting	Alum Solution

**Table XVII (continued)      Chemicals which are compatible  
with GRAFOIL Flexible Graphite**

**Chemical Reagents**

Aluminum Chloride	Asphalt
Aluminum Fluoride	Barium Carbonate
Aluminum Hydroxide (Boehmite)	Barium Chloride (Aqueous)
Aluminum Sulfate	Barium Hydroxide
Aluminum Sulphate, Aqueous	Barium Nitrate (Aqueous)
Amine	Barium Sulfate
Amine, Fat Condensate	Barium Sulfide
Ammonia	Beer
Ammonia Liquid	Beer Wort
Ammonia, Anhydrous	Beet Juice
Ammonia, Aqua	Beet Pulp
Ammonia, Aqueous	Beet Sugar Solution
Ammonia, Gas (Dry)	Benoic Acid Sol
Ammonium Bicarbonate	Benzaldehyde
Ammonium Bifluoride	Benzene
Ammonium Bisulfate	Benzene Hexachloride
Ammonium Carbonate	Benzoic Acid
Ammonium Chloride	Benzol
Ammonium Hydroxide	Benzyl Sulfonic Acid
Ammonium Nitrate	Black Liquor
Ammonium Persulfate	Bleach Solution (0-25% conc. @R.T.)
Ammonium Phosphate	Boiler Feed Water
Ammonium Sulfate	Borax (Sodium Borate)
Ammonium Thiocyanate	Borax Solution
Amyl Acetate	Boric Acid
Amyl Alcohol	Boric Acid, Aqueous
Amyl Nitrate	Boron Trichloride
Aniline	Brine Calcium
Aniline Hydrochloride	Brine, Calcium & Sodium Chloride
Aniline Hydrochlorine	Brine, Chloride PH8
Anti-Foaming Agent	Brine & Magnesium Chloride
Aromatic Fuels	Brine, Sea Water
Aniline Hydrochloride	Bromo Methane
Arsenic Acid	Bunker C Fuel Oil
Arsenic Trichloride	Butadiene

**Table XVII (continued)      Chemicals which are compatible  
with GRAFOIL Flexible Graphite**

**Chemical Reagents**

Butane	Caustic Zine Chloride
Butane (LPG)	Chloracetone
Butanol	Chloral Hydrate
Butyl Acetate	Chloride of Lime
Butyric Acid	Chloride of Zinc
Butyl Alcohol	Chlorinated Solvents
Butylamine	Chlorine
Butyl Caritol	Chlorine (Dry and Wet)
Butylene	Chlorine Water
Butylene (Butene) (Ethylethylene)	Chloroacetic Acid
Butyl Phthalate	Chlorobenzene
Calcium Carbonate	Chloroethylbenzene
Calcium Chloride	Chloroform
Calcium Hydroxide	Chloropicrin
Calcium Hypochlorite	Chlorosulfonic Acid
Calcium Magnesium Chloride	Citric Acid
Calcium Phosphate	Citric Acid, Aqueous
Calcium Sulfate	Cocoa Butter
Cane Juice	Coconut Acid, Fatty
Carbolic Acid	Coconut Oil
Carbonate of Soda	Cod Liver Oil
Carbon Bisulfide	Condensate (Water)
Carbon Dioxide	Cooling Water
Carbon Disulfide	Copper Acetate
Carbonic Acid	Copper Ammonium Acetate
Carbonic Acid, Aqueous	Copperas
Carbon Monoxide	Copper Cyanide
Carbon Tetrachloride	Copper Nitrate
Carbon Tetrachloride (Anhydrous)	Copper Sulfate
Catsup	Corn Oil
Caustic	Creosol, META
Caustic Cyanogen	Cresols
Caustic Potash	Creosote (Coal Tar)
Caustic Soda	Cresylic Acid
Caustic Strontia	Crude Oil
Caustic Sulfide	Cupric Chloride

**Table XVII (continued)      Chemicals which are compatible  
with GRAFOIL Flexible Graphite**

**Chemical Reagents**

Cupric Sulfate	Dinitrochlorobenzene & Styrene
Cupros Ammonia Acetate	Dinitrotoluene
Cutting Oil	Diethyl Phthalate
Cyanogen in Water	Diethyl-Amine
Cyclohexene	Diphenyl
Cyclohexanone	Dow Corning Silicone Fluid
DDT Solution (Kerosine Solv)	Dowtherm (All Types)
DDT Solution (Toluene Solv)	Dye Wood Liquor
De-Butanizer Reflux	Esso Therm 500
De-Ethanizer Charge	Ethane
Demineralized Water	Ethanol (Ethyl Alcohol)
Denatured Alcohol	Ethanolamine
De-Propanizer Reflux	Ether
Diacetone Alcohol	Ethyl Acetate
Dibromoethyl Benzene	Ethyl Acrylate
Dibutyl Amine	Ethyl Alcohol (Ethanol)
Dibutyl Cellusolve Adipate	Ethyl Benzoate
Dibutyl Phthalate	Ethyl Cellosolve
Dibutylether	Ethyl Chloride
Dichlorobenzene	Ethyl Chlorocarbonate
Dichloroethane	Ethyl Chloroformate
Diesel Fuel	Ethyl Chlorohydrin
Diethanolamine (Dea)	Ethylene (Ethene)
Diethyl Carbonate	Ethylene Chloride
Diethyl Ether	Ethylene Chlorohydrin
Diethylamine	Ethylene Diamine
Diethylene Glycol	Ethylene Bromide
Diethylene Triamine	Ethylene Dichloride
Di-Isobutyl Ketone	Ethylene Glycol (anti-freeze)
Di-Isopropyl Ketone	Ethylene Oxide
Dimethyl Formaldehyde	Ethylene Oxide & Freon 12
Dimethyl Hydrazene	Ethylene Trichloride
Dimethyl Phthalate	Ethyl Ether (Ethyl Oxide)
Dimethyl Terephthalate	Ethyl Formate
Dinitrochlorobenzene	Ethyl Hexanol

**Table XVII (continued) Chemicals which are compatible with GRAFOIL Flexible Graphite**

**Chemical Reagents**

Ethyl Mercaptain-water	Gasoline Aromatic
Ethyl Oxide	Gasoline Hi-Test w/Mercaptan, H <sub>2</sub> S
Ethyl Pyridine	Gelatin
Fatty Acids	Glaubers Salt
Fatty Acids, Oleic	Glycerine (Glycerol)
Ferric Chloride	Glycols
Ferric Chloride (Aqueous)	Grease
Ferric Nitrate	Glue Sizing
Ferric Sulfate (Aqueous)	Green Sulphate Liquor
Ferrous Chloride	Halon
Ferrous Sulfate	Heat Transfer Fluids (Petroleum-Oil Based)
Flourolube	Helium
Fluoboric Acid	Heptane
Fluorine	Hexachloro Acetone
Fluosilicic Acid	Hexane
Folic Acid	Hexanol
Formaldehyde (Methanol)	Hexene (Butylethylene)
Formalin	Hexone
Formic Acid	Hexyl Alcohol (Hexanol)
Freon 11 & Refrig. Oil	Hops
Freon 113 & Refrig. Oil	Hydrazine
Freon 114 & Refrig. Oil	Hydrobromic Acid
Freon 12 & Refrig. Oil	Hydrochloric Acid
Freon 121 & Refrig. Oil	Hydrocyanic Acid
Freon 22 & Refrig. Oil	Hydroflouric Acid
Freons, Liquid	Hydrogen
Fruit Acid	Hydrogen Chloride
Fruit Juices	Hydrogen Fluoride
Fuel Oil	Hydrogen Sulfide
Fuel Oil #6	Hydrogen Sulfide-water
Fuel Oil, Acidic	Hypochlorous Acid
Furfural	Hydrogen Chloride
Furnace Oil	Insecticides (Aromatic)
Gas Oil	Insecticides (Nonaromatic)
Gasoline	Iodine
Gasoline 100 & 130 Oct.	Iodoform

**Table XVII (continued)      Chemicals which are compatible  
with GRAFOIL Flexible Graphite**

**Chemical Reagents**

Iso-Octane	Lithium Chloride (Aqueous)
Isobutane	Lye, Caustic
Isobutyl Alcohol	Lye, Salty
Isobutyl Methyl Ketone	Magnesium Chloride (Bischofite)
Isobutylene	Magnesium Hydroxide (Brucity)
Isodecane	Magnesium Sulfate (Epsom Salt)
Isopentane	Magnesium Sulfate (Aqueous)
Isopropanol	Magnesium Sulfite (Aqueous)
Isopropyl Acetate	Maleic Anhydride
Isopropyl Alcohol	Maleic Hydrazide
Isopropyl Chloride	Manganese Chloride (Aqueous)
Isopropyl Ether	Manganese Sulfate (Aqueous)
Isopropylamine	Marsh Gas (Methane)
Isopropyl Acetone	Marsh, Anti-Biotic Fermentation, No Solvent
JP-3	Marsh, With Solvent
JP-4	Mayonnaise
JP-5	MEA (Monoethanolamine)
JP-6	MEA- With Copper Sulfate
JP-X	Melamine Resins
Kerosene	Mercaptan
Ketones	Mercuric Chloride
Lacquer (MEK Sol)	Mercury
Lacquer Thinners	Mercury Salts
Lactic Acid	Mercury Vapors
Lard (Animal Fat)	Mesityl Oxide (Ketone)
Latex	Methane (Marsh Gas)
Lavender Oil	Methanol (Methyl Alcohol)
Lead Acetate (Liquid)	Methyl Acetate
Lead Nitrate	Methyl Acrylate
Lead Sulphamate	Methyl Alcohol (Methanol)
Lindol	Methyl Bromide (Bromomethane)
Linseed Oil	Methyl Butyl Ketone
Liquor, Steep	Methyl Cellosolve
Liquor, Sulphate	Methyl Chloride (Chloromethane)
Lithium Bromide Brine	Methyl Chloride (Anhydrous)

**Table XVII (continued)      Chemicals which are compatible  
with GRAFOIL Flexible Graphite**

**Chemical Reagents**

Methyl Cyclopentane	Nitrochloroform
Methyl Dichloride	Nitrobenzene
Methyl Ether	Nitrobenzine
Methyl Ethyl Ketone (MEK)	Nitroethane
Methyl Formate	Nitromethane
Methyl Isobutyl Ketone	Nitrogen Gas
Methyl Isopropyl Ketone	Nitropropane
Methylene Chloride (Dichloromethane)	Oakite
Methylene Dichloride	Oil & Ammonia
Milk	Oil, Animal, Bone
Milk of Lime	Oil, Animal, Cod
MIL F-25558 (RJ-1)	Oil, Animal, Lard
MIL H-5606 (HFA)	Oil, Animal, Menhadden
MIL H-5606 (J43)	Oil, Animal, Neatsfoot
MIL L-7808	Oil, Animal, Sperm
MIL O-8200 (Hydr)	Oil, Animal, Whale
MIL O-8515	Oil, Bunker 'C'
Mine Water	Oil, Coal Tar
Mineral Spirits	Oil, Creosote, Sweet
Miscella 20% Soya Oil	Oil, Crude, Sweet
Mobiltherm	Oil, Diesel, #2D
Molasses	Oil, Diesel, #3D
Monochlorobenzene	Oil, Diesel, #4D
Monoethanolamine (MEA)	Oil, Diesel, #5D
Naphtha	Oil, Essential
Naphtha Crude	Oil, Fed. Spec #9170, #9250
Napthalene	#9370, #9500
Natural Gas Liquid	Oil, Fed. Spec #10, #20, #30
Neatsfoot Oil	Oil, Fed. Spec SAE 20, 30, 40, 50, 60,
Nickel Acetate	70, 90, 140, 250
Nickel Chloride	Oil, Fuel #1, #2, #3, #5A, #5B, #6
Nickel Cobalt Sulfate, 5% H <sub>2</sub> SO <sub>4</sub>	Oil, Insulating
Nickel Salts	Oil, Kerosene
Nickle Sulfate	Oil, Lean
Nicotine Sulfate	Oil, Linseed (Raw)



**Table XVII (continued)      Chemicals which are compatible  
with GRAFOIL Flexible Graphite**

**Chemical Reagents**

Oil, Lubricating #8	Petroleum Ether
Oil, Lubricating Diesel #9110	Phenol (Carbolic Acid)
Oil, Mineral Lard Cutting, Fed. Spec. #1	Phenol, Formaldehyde Mix
Oil, Mineral Lard Cutting, Fed. Spec. #2	Phenyl Acetic Acid
Oil, Navy Spec., Navy II	Phidolene
Oil, Quenching	Phosphoric Acid
Oil, Rich	Phosphorous Trichloride
Oil, Turbine Lube	Photographic Developers
Oil, Vegetable, Castor	Phthalic Anhydride
Oil, Vegetable, China Wood	Phthalic Esters
Oil, Vegetable, Coconut	Picric Acid, Molten
Oil, Vegetable, Corn	Picric Acid, Water Solution
Oil, Vegetable, Cottonseed	Plasticizer
Oil, Vegetable, Linseed (Raw)	Plating Solutions (Not Chrome)
Oil, Vegetable, Olive	Poly Glycols
Oil, Vegetable, Palm	Poly Vinyl Acetate
Oil, Vegetable, Peanut	Potable Water
Oil, Vegetable, Rape Seed	Potash (Plant Liquor)
Oil, Vegetable, Rosin	Potash Alum
Oil, Vegetable, Sesame	Potash, Sulfide
Oil, Vegetable, Soya Bean	Potassium Bicarbonate
Olefin, Crude	Potassium Bromide
Oleums	Potassium Carbonate
Orthodichloro Benzene	Potassium Chloride
Oxygen, Liquid and Gaseous (below 600°F)	Potassium Cyanides
Paracymene	Potassium Hydroxide
Paradichlorobenzene	Potassium Perfluoro Acetate
Paraffin, Liquid	Potassium Permanganate
Pectin, Liquor	Potassium Phosphate
Penicillin, Liquid	Potassium Sulfate
Pentachlorophenol	Potassium, Bicarbonate
Pentane	Propane
Pentasol	Propiolactone, Beta
Perchloroethylene	Propionaldehyde
Petrolatum	Propyl Alcohol

**Table XVII (continued)      Chemicals which are compatible  
with GRAFOIL Flexible Graphite**

**Chemical Reagents**

Propylene (Propane)	Sodium Phosphate, Tribasic
Propylene Glycol	Sodium Plumbite
Propylene Oxide	Sodium Silicate
Pulp Stock	Sodium Sulfate
Pyridine	Sodium Sulfide
Pyrogallic Acid	Sodium Tetraborate
Pyroligneous Acid	Sodium Thiosulfate
Raffinate	Solvasol 1, 2, 3, 73, 74
Refrigerants 11, 12, and 134A	Sorbitol
Sal Ammoniac (Ammonium Chloride)	Stannic Chloride
Salt Cake	Starch
Salt Water	Steam
Sea Water	Steam, Superheated
Sewage	Stearic Acid & Oleic
Shellac	Steric Acid (Octadecandic Acid)
Silicone Oils and Greases	Stoddard Solvent
Silver Nitrate	Strontium Nitrate
Skydrol 500, 7000	Styrene (Monomer) (Vinylbenzene)
Soap, Liquors	Sulfate of Lime
Soap, Solutions	Sulfide of Hydrogen
Sodium Acetate	Sulfide of Sodium
Sodium Aluminate	Sulfite Pulp
Sodium Bicarbonate	Sulfur Chloride
Sodium Bisulfate	Sulfur Dioxide
Sodium Bisulfite	Sulfur Monochloride
Sodium Borate	Sulphonated Fatty Alcohol
Sodium Chloride	Sulphonated Vegetable Oils
Sodium Chloride Sol ( Salt)	Sulphuric Chlorohydrin
Sodium Cyanamide	Syrup
Sodium Hydroxide	Tall Oil
Sodium Hydrosulfite	Tallow
Sodium Metasilicate	Tannic Acid
Sodium Perborate	Tanning Liquors
Sodium Phosphate, Dibasic	Tar & Ammonia w/Water
Sodium Phosphate, Meta	Tar, Bituminous
Sodium Phosphate, Mono	

**Table XVII (continued)      Chemicals which are compatible  
with GRAFOIL Flexible Graphite**

**Chemical Reagents**

Tar, Pine	Water, Borated
Tartaric Acid, Aqueous	Water, Brackish
Terephthalic Acid	Water, Clean Untreated
Tetrachloroethane	Water, Condensate
Tetrachloroethylene	Water, Cooling Tower
Tetrahydrofuran	Water, Deaerated
Tetraphenyl	Water, Distilled
Therminol (All Types)	Water, Fresh
Therminol #1, 2 & 3	Water, Heavy
Titanium Tetrachloride	Water, Hot
Toluene (Toluol) (Methylbenzene)	Water, Mine
Tomato Pulp	Water, River
Toxaphene	Water, s/Sol. Oil
Trichlorobenzene	Water, Salt & Sea, Solution
Trichloroethane	Water, Soapy
Trichloroethylene	Whiskey
Trichloronitromethane	White Liquor
Tricresylphosphate	White Water, Paper Mill
Triethylamine	Wine
Trifluorovinylchloride	Wood Pulp (Stock)
Trisodium Phosphate	Wood Vinegar
Turpentine	Wort (Beer Wort)
Ucon (All types)	Xylene (Dimethylbenzene)
Urea	Yeast
Varnish, Aromatic	Zeolite Treated Water
Varnish, Non-Aromatic	Zinc Ammonium Chloride
Vegetable Juices	Zinc Chloride
Vetrocoke Solution	Zinc Cyanide
Vinegar	Zinc Nitrate
Vinyl Chloride	Zinc Phosphate
Vinyl Pyridine	Zinc Sulfate
Vinylidene Chloride	
Water, Boiler Feed	