



# High Temperature 4570 Vessel



Parr Instrument Company has designed a vessel based on the 4570 model that is rated 5000 psi at 550 °C. This high temperature vessel maintains the 3.75" ID from the 4571/4572 product family and features metal-to-metal or metal-gasketed connections for all head fittings to ensure leak tight operation at higher temperature conditions.

The vessel illustrated to the left has a 1 Liter internal volume and is customized for the higher temperature rating with special split ring and bolts. A custom head configuration can be requested based on customer preference.

**Let Us Build One For You.**

**Parr Instrument Company** • 211 Fifty-Third Street, Moline, Illinois 61265-1770 U.S.A.

Phone: 1-800-872-7720 or (309) 762-7716 • Fax: (309) 762-9453 • [parr@parrinst.com](mailto:parr@parrinst.com) • [www.parrinst.com](http://www.parrinst.com)